


# Product Information Sheet



<b>Item Name</b>	<b>Breast Tumor Cell Line Pellet, RNA Preserved</b>	<b>Prod. Code:</b>	<b>2250-0055</b>
<b>Stage / Grade</b>	Cell Pellet	<b>Optional Description:</b>	ZR-75-1 Cell Pellet
<b>Source:</b>	Human Cell Line-ZR-75-1		

	<u>Method Used:</u>	<u>Acceptable Result:</u>
<b>Protein Quantitation :</b>	N/A	N/A mg
<b>Product Quality:</b>	N/A	N/A
<b>Total Protein:</b>	N/A	N/A

**Buffer Matrix & pH:** N/A

**Filtration:** N/A

**Optimum Storage Conditions:** Short term : -20°C Long term: -70°C

**Shelf Life:** TBD (stored at recommended long term temperature)

**Other:** N/A

## Hazard Testing & Warnings:

<b>Viral Testing:</b>	The cell line supernatant tested negative for HIV1,2, HIVag, HbsAg, and HCV. Please note that the samples tested were cell culture supernatants, not serum, and therefore not appropriate for use in these assays.
<b>Caution Statement:</b>	<b>NO TEST CAN GUARANTEE A PRODUCT TO BE NON-INFECTIOUS. THEREFORE, ALL MATERIAL DERIVED FROM HUMAN FLUIDS OR TISSUES SHOULD BE CONSIDERED AS POTENTIALLY INFECTIOUS.</b>
<b>Use Statement:</b>	<b>-For use in research or remanufacture of non-injectable products not subject to license under PHS Act Sec. 351 only.</b>