

Product Information Sheet



Item Name	Ovar. Cancer Antigen, PP, LCR	Prod. Code:	8375-1303
Stage / Grade	Partially pure, LCR	Optional Description:	Ovarian Cancer Antigen, Partially Purified, Low Cross Reactivity
Source:	Human Fluids		

	Method Used:	Acceptable Result:
Protein Quantitation :	Centocor Ca-125 Assay	50,000 - 200,000 U/ml U/ml
Product Quality:	U/ml ÷ OD280 =relative purity.	≥5000
	Centocor CA 15-3 RIA Assay	<25% Ca-125 result
	Centocor 19-9 RIA Assay	<25% Ca-125 result
Total Protein:	OD280	≤20 OD units

Buffer Matrix & pH: PBS Buffer with 0.1% NaN₃, pH=7.4±0.2

Filtration: 0.2u

Optimum Storage Conditions: Short term : 4°C Long term: -20°C

Shelf Life: 5 years at -20°C (stored at recommended long term temperature)

Other: N/A

Hazard Testing & Warnings:

Viral Testing: All donor units tested negative for HIV1,2, HIVag, HbsAg, and HCV.

Caution Statement: NO TEST CAN GUARANTEE A PRODUCT TO BE NON-INFECTIOUS. THEREFORE, ALL MATERIAL DERIVED FROM HUMAN FLUIDS OR TISSUES SHOULD BE CONSIDERED AS POTENTIALLY INFECTIOUS.

Use Statement: -For use in research or remanufacture of non-injectable products not subject to license under PHS Act Sec. 351 only.