

# Product Information Sheet



<b>Item Name</b>	<b>RF Control Serum</b>	<b>Prod. Code:</b>	<b>8875-3000</b>
<b>Stage / Grade</b>	Control	<b>Optional Description:</b>	Rheumatoid Factor Control Serum
<b>Source:</b>	Human serum		

	<u>Method Used:</u>	<u>Acceptable Result:</u>
<b>Protein Quantitation :</b>	Behring RF Latex or Behring RF Nephelometric Test	1,100-1,200 IU/ml IU/ml
<b>Product Quality:</b>	OD 570 and OD 660	OD 570 0.06 , OD 660 0.07
<b>Total Protein:</b>	OD280	≤ 40 OD Units

**Buffer Matrix & pH:** 0.1M Tris, 0.2M NaCl, pH 7.5, with 0.1% NaN3

**Filtration:** 0.2u

**Optimum Storage Conditions:** Short term : 4°C Long term: -20°C

**Shelf Life:** 1 year at -20°C (stored at recommended long term temperature)

**Other:** n/a

## Hazard Testing & Warnings:

**Viral Testing:** All donor units tested negative for HIV1,2, HIVag, HbsAg, and HCV.

**Caution Statement:** NO TEST CAN GUARANTEE A PRODUCT TO BE NON-INFECTIOUS. THEREFORE, ALL MATERIAL DERIVED FROM HUMAN FLUIDS OR TISSUES SHOULD BE CONSIDERED AS POTENTIALLY INFECTIOUS.

**Use Statement:** -For Research or Remanufacturing Use Only